



**BARBOUR**  
PLASTIC SURGERY

**JP Drain Record**

Drain #1			Drain #2		
Date		Output (cc/ml)	Date		Output (cc/ml)
	am			am	
	pm			pm	
	total			total	
	am			am	
	pm			pm	
	total			total	
	am			am	
	pm			pm	
	total			total	
	am			am	
	pm			pm	
	total			total	
	am			am	
	pm			pm	
	total			total	
	am			am	
	pm			pm	
	total			total	

### **Drain Record Information/Instructions:**

- Drains are placed to prevent excess fluid from accumulating at the surgery site.
- The drain must be emptied 2-3 times a day or when 2/3 full.
- Pour out the drainage into a clean measuring cup, record the amount of drainage (in mL or CCs) and the time of day and dispose of drainage in the toilet and flush.
- After the fluid is emptied, squeeze the bulb tight and recap to maintain the suction. The bulb suction is working when the bulb is flat or deflated. If blockages arise in the drain hold tubing securely and pinch and slide fingers down tubing to release.
- If you have more than one drain in place, please indicate which drain is being recorded on your record sheet.
- It is ok to shower while drains are in place, being very careful to not let drains hang down or pull at drain sites. Do not soak in the tub or submerge the drain sites under water. Wash gently around drain sites with gentle soap and water.
- Check skin around drain, slight redness is normal. If there is increasing redness, pain or leakage around drain, contact Dr. Barbour's office.